

Product Datasheet

LAMP2 antibody orb33407

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LAMP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44961
Clone ID	RB13276
Uniprot ID	P13473
NCBI	3920, NP_002285.1, NP_001116078.1, NP_054701.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

