




Product Datasheet

LC3C (phospho-Ser12) antibody orb34558

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	MAP1LC3A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14272
Clone ID	RB10838, RB16769
Uniprot ID	Q9H492
NCBI	84557 , NP_852610.1 , NP_115903.1
Dilution Range	WB: 1:1000, DOT: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of CHO cell culture using LC3C (phospho-Ser12) antibody (primary antibody at 1:1000)



Western blot analysis of level of LC3 using LC3C (phospho-Ser12) antibody (dilution of primary antibody at 1:1000)