

Product Datasheet

LC3 (phospho-Thr12) antibody orb34767

Description	Rabbit polyclonal antibody to LC3 (phospho-Thr12)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MAP1LC3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14688
Clone ID	RB12679
Uniprot ID	Q9GZQ8
NCBI	81631 , NP_073729.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

