

Product Datasheet

Mouse Anti-Human IgE antibody orb349365

Description	Mouse monoclonal antibody to Mouse Anti-Human IgE.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only.
Isotype	IgG1,?
Clonality	Monoclonal
Clone ID	1497CT272.23.66
Dilution Range	ELISA: 1:5000
Expiration Date	12 months from date of receipt.

