

## Product Datasheet

### Mouse Anti-Human IgE antibody orb349366

<b>Description</b>	Mouse monoclonal antibody to Mouse Anti-Human IgE.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Form/Appearance</b>	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG1,?
<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1497CT744.79.17
<b>Dilution Range</b>	ELISA: 1:5000
<b>Expiration Date</b>	12 months from date of receipt.

