

Product Datasheet

APG3 (phosphorylated) antibody orb34972

Description	Rabbit polyclonal to APG3 (phosphorylated)
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ATG3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35864
Clone ID	RB18431
Uniprot ID	Q9NT62
NCBI	64422 , NP_001265641.1 , NP_071933.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

