



Product Datasheet

BTG1 antibody orb35408

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	BTG1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19209
Clone ID	RB26174
Uniprot ID	P62324
NCBI	694 , NP_001722.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 cell line lysates (35ug/lane) using BTG1 antibody (primary antibody dilution at: 1:1000)