



Product Datasheet

SIRT7 antibody orb36407

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SIRT7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44898
Clone ID	RB1973
Uniprot ID	Q9NRC8
NCBI	51547 , NP_057622.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate using SIRT7 antibody (primary antibody dilution at: 1:1000)