



Product Datasheet

Annexin V antibody orb36729

Description	Rabbit polyclonal antibody to Annexin V
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ANXA5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35937
Clone ID	RB19462
Uniprot ID	P08758
NCBI	308 , NP_001145.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin-fixed paraffin embedded human hepatocarcinoma tissue using Annexin V antibody