




Product Datasheet

CD163 antibody orb37792

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CD163
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	125451
Clone ID	RB18866
Uniprot ID	Q86VB7
NCBI	9332 , NP_004235.4 , NP_981961.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of CD163 antibody in K562 cell line lysates



Immunocytochemistry on J774 murine macrophages (DAPI image).