



Product Datasheet

IGFBP3 antibody orb38296

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse
Predicted Reactivity	Bovine, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	This IGFBP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 162-189 amino acids from human IGFBP3.
Target	IGFBP3
Form/Appearance	PBS with 0.09% (W/V) sodium azide
Storage	Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Peptide Affinity purified
MW	31.674 kDa
Uniprot ID	P17936 (Human)
Entrez	3486
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of
CHO
tissue lysates
(15ug/lane) using IGFBP3
antibody (primary antibody
dilution at: 1:1000)