



Product Datasheet

Angiotensinogen antibody orb38626

Description	Rabbit polyclonal antibody to Angiotensinogen
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	AGT
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53154
Clone ID	RB17036
Uniprot ID	P01019
NCBI	183, NP_000020.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in HepG2 cell line lysates (35ug/lane) using Angiotensin II antibody (primary antibody dilution at: 1:1000)