




Product Datasheet

CAV2 antibody orb39333

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CAV2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18291
Clone ID	RB16483
Uniprot ID	P51636
NCBI	858, NP_937855.1, NP_001193676.1, NP_001224.1, NP_001193677.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MDA-MB231 cell using CAV2 antibody (primary antibody dilution at: 1:10-50)



Immunohistochemical staining of paraffin embedded human skeletal muscle tissue using CAV2 antibody (primary antibody dilution at: 1:50-100)