

Product Datasheet

DKK1 antibody orb40298

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DKK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28672
Clone ID	RB21368
Uniprot ID	O94907
NCBI	22943 , NP_036374.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

