

Product Datasheet

PSMG2 antibody orb40398

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Immunogen is PSMG2 (C-term).
Target	PSMG2
Form/Appearance	Supplied in PBS with 0.09% (W/V) sodium azide.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Peptide affinity purified.
Uniprot ID	Q969U7
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

