



Product Datasheet

HDAC2 antibody orb40465

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	HDAC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55364
Clone ID	RB24526
Uniprot ID	Q92769
NCBI	3066 , NP_001518.3
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cells using HDAC2 antibody (primary antibody dilution at: 1:10-50)



Immunohistochemical staining of paraffin embedded human lung carcinoma tissue using HDAC2 antibody (primary antibody dilution at: 1:50-100)