




Product Datasheet

Gastrin Receptor antibody orb40494

Description	Rabbit polyclonal antibody to Gastrin Receptor
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GRPR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43199
Clone ID	RB24634
Uniprot ID	P30550
NCBI	2925, NP_005305.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin fixed and paraffin embedded brain tissue using GRPR antibody



Flow cytometric analysis of GRPR antibody in HL-60 cell line lysates