

Product Datasheet

FGFR1-Y307 antibody orb40518

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Gallus, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	This FGFR1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 285-314 amino acids from human FGFR1.
Target	FGFR1-Y307
Form/Appearance	PBS with 0.09% (W/V) sodium azide
Storage	Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Peptide Affinity purified
MW	91.868 kDa
Uniprot ID	P11362 (Human)
Entrez	2260
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

