




Product Datasheet

IMP3 antibody orb40603

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IMP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21850
Clone ID	RB24885
Uniprot ID	Q9NV31
NCBI	55272 , NP_060755.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin-fixed paraffin embedded mouse testis tissue using IMP3 antibody



Flow cytometric analysis of HeLa cells (right histogram) compared to a negative control cell (left histogram) using IMP3 antibody