




Product Datasheet

ID4 antibody orb40682

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ID4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16622
Clone ID	RB24654
Uniprot ID	P47928
NCBI	3400, NP_001537.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin embedded human bladder carcinoma tissue using ID4 antibody (primary antibody dilution at: 1:50-100)



Western blot analysis of in mouse bladder tissue lysates (35ug/lane) using ID4 antibody (primary antibody dilution at: 1:1000)