

## Product Datasheet

### GFAP antibody orb97669

<b>Description</b>	Mouse monoclonal to GFAP. This GFAP monoclonal antibody is generated from mouse immunized with GFAP recombinant protein.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, IHC-P, WB
<b>Target</b>	GFAP. This GFAP monoclonal antibody is generated from mouse immunized with GFAP recombinant protein.
<b>Form/Appearance</b>	Supplied in PBS with 0.09% (W/V) sodium azide.
<b>Storage</b>	Maintain refrigerated at 2-8°C. for up to 6 months. For long term storage store at -20°C. in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG2b,k
<b>Clonality</b>	Monoclonal
<b>Purity</b>	This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
<b>MW</b>	49880 Da
<b>Clone ID</b>	183CT3.1.5
<b>NCBI</b>	<a href="#">NP_001124491.1</a>
<b>Dilution Range</b>	WB: 1:100-500, IHC-P: 1:50-100, IF/ICC: 1:10-50
<b>Expiration Date</b>	12 months from date of receipt.

