

Product Datasheet

TGFB2 antibody orb97672

| | |
|----------------------------|--|
| Description | Mouse monoclonal to TGFB2. This TGFB2 monoclonal antibody is generated from mouse immunized with TGFB2 recombinant protein. |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, ICC, IF, IHC-P, WB |
| Immunogen | This TGFB2 monoclonal antibody is generated from mouse immunized with TGFB2 recombinant protein. |
| Target | TGFB2. This TGFB2 monoclonal antibody is generated from mouse immunized with TGFB2 recombinant protein. |
| Form/Appearance | PBS with 0.09% (W/V) sodium azide |
| Storage | Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only. |
| Isotype | IgG1 Kappa |
| Clonality | Monoclonal |
| Purity | This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS. |
| MW | 47.748 kDa |
| Clone ID | 220ct16.4.3.1 |
| Uniprot ID | P61812 (Human) |
| NCBI | NP_001129071.1 |
| Entrez | 7042 |
| Dilution Range | WB: 1:100-1000, IHC-P: 1:10-50, IF/ICC: 1:10-50 |
| Expiration Date | 12 months from date of receipt. |

