

Product Datasheet

STAT3 antibody orb97688

Description	Rabbit polyclonal to STAT3. This STAT3 antibody is generated from rabbit immunized with a KLH conjugated synthetic peptide from the C-terminal region of human STAT3, aa 711-743.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	STAT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88068
Clone ID	RB38197
Uniprot ID	P40763
NCBI	6774 , NP_998827.1 , NP_003141.2 , NP_644805.1
Dilution Range	WB: 1:1000, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.

