

Product Datasheet

SLCO1B1 antibody orb613766

Description	Goat polyclonal antibody to SLCO1B1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, P-ELISA
Immunogen	C-KYQEKDINASEN
Target	SLCO1B1
Form/Appearance	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Antigen affinity chromatography using the immunizing peptide.
MW	76.4
NCBI	NP_006437.3
Entrez	10599
Dilution Range	P-ELISA: 1:32000,IF: 10 µg/ml
Expiration Date	12 months from date of receipt.
Application Notes	Notes: Immunofluorescence: Strong expression of the protein seen in the plasma membrane of HepG2 cells. Recommended concentration: 10µg/ml. This antibody has been successfully used in IF on Human: Malinen MM et al (2019). PMID: 30685239.