


Product Datasheet

Rat TROY ELISA Kit orb397145

Description	Sandwich ELISA for quantitative detection of Rat TROY concentrations.
Reactivity	Rat
Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	30 pg/ml
Tested Applications	ELISA
Target	TNFRSF19
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	D3ZUJ0
Entrez	290300
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates



Representative Standard Curve