

Product Datasheet

Human C1QTNF6 ELISA Kit orb404202

Description	Human Complement C1q tumor necrosis factor-related protein 6(C1QTNF6) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Target	C1QTNF6
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BXI9
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm