

Product Datasheet

Affinipure Rabbit anti Goat IgG(H+L) antibody (Gold) orb633274

Description	Rabbit polyclonal antibody against Affinipure Rabbit anti Goat IgG(H+L) conjugated to Gold
Conjugation	Gold
Form/Appearance	liquid
Preservatives	1%rAlbumin
Storage	2mM boric acid,1%trehalose, 0.1%NaN3, pH9.0, 4? for 24 months
Note	For research use only.
Expiration Date	12 months from date of receipt.