

## Product Datasheet

### Rabbit anti Human IgM antibody (HRP) orb633275

|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody against Rabbit anti Human IgM conjugated to HRP                                                          |
| <b>Conjugation</b>         | HRP                                                                                                                                 |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                           |
| <b>Form/Appearance</b>     | liquid                                                                                                                              |
| <b>Preservatives</b>       | 3mg/ml rAlbumin ( IgG-Free, Protease-Free)                                                                                          |
| <b>Storage</b>             | 0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20° for 24 months , avoid repeated freezing and thawing. |
| <b>Note</b>                | For research use only.                                                                                                              |
| <b>Dilution Range</b>      | 1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                     |