

Product Datasheet

Rabbit anti Human IgG antibody (HRP) orb633277

Description	Rabbit polyclonal antibody against Rabbit anti Human IgG conjugated to HRP
Conjugation	HRP
Tested Applications	ELISA, WB
Form/Appearance	liquid
Concentration	3mg/ml BSA (IgG-Free, Protease-Free)
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20° for 24 months , avoid repeated freezing and thawing.
Note	For research use only.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.