

Product Datasheet

Human ANPEP ELISA Kit orb638773

Description	Human alanyl (membrane) aminopeptidase (ANPEP)ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Sensitivity	7.81 pg/mL
Tested Applications	ELISA
Target	ANPEP
Storage	Store at 2-8°C. Please refer to protocol.
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P15144
Expiration Date	Please enquire.
Application Notes	Sample Types: serum, plasma, cell culture supernates Assay Time: 1-5h Sample Volume: 50-100ul Detection Wavelength: 450 nm