

Product Datasheet

C5b-9 antibody {CF647} orb100655-CF647

Description	Rabbit polyclonal antibody to C5b-9
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human C5b-9. Please contact us for the exact immunogen sequence. The peptide is available as orb12903 .
Target	C5b-9
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 520 µl): 0.192 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	188 kDa
Uniprot ID	P01031 (Human)
Dilution Range	IF/ICC: 1:100-800
Expiration Date	12 months from date of receipt.