

Product Datasheet

CCDC25 antibody {CF647} orb2517-CF647

Description	Rabbit polyclonal antibody to CCDC25
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Equine, Porcine, Sheep
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa81-160 of human CCDC25
Target	CCDC25
Form/Appearance	Liquid: PBS, 0.2 mM sodium Azide
Concentration	- 100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	24 kDa
Uniprot ID	Q86WR0 (Human); Q78PG9 (Mouse)
Entrez	55246
Dilution Range	IF: 1:100-500
Expiration Date	12 months from date of receipt.