

Product Datasheet

Cathepsin G antibody {CF647} orb500931-CF647

Description	Rabbit polyclonal antibody to Cathepsin G
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa51-110 of human Cathepsin G
Target	Cathepsin G
Form/Appearance	Liquid: PBS
Concentration	- 100 µg (in 522 µl): 0.19 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	28 kDa
Uniprot ID	P08311 (Human)
Entrez	1511
Dilution Range	IF: 1:100-500
Expiration Date	12 months from date of receipt.