

Product Datasheet

Collagen IV antibody {CF405M} orb6174-CF405M

Description	Rabbit polyclonal antibody to Collagen IV
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Equine, Porcine
Conjugation	CF405M
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived between 1411-1490 amino acids of human Collagen IV
Target	Collagen IV
Form/Appearance	Liquid: PBS
Concentration	- 100 µg (in 410 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	165 kDa
Uniprot ID	P02462 (Human)
Entrez	1282
Dilution Range	IF: 1:100-500
Expiration Date	12 months from date of receipt.