

Product Datasheet

CROCC antibody {CF647} orb182754-CF647

Description	Rabbit polyclonal antibody to CROCC
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Equine, Porcine
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa1861-2007 of human CROCC
Target	CROCC
Form/Appearance	Liquid: PBS
Concentration	- 100 µg (in 410 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	228 kDa
Uniprot ID	Q5TZA2 (Human)
Entrez	9696
Dilution Range	ICC: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.