

Product Datasheet

Neutrophil Elastase antibody {CF647} orb500877-CF647

Description	Rabbit polyclonal antibody to Neutrophil Elastase
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	CF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa61-140 of human Neutrophil Elastase
Target	Neutrophil Elastase
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	- 100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	26 kDa
Uniprot ID	P08246 (Human); Q3UP87 (Mouse)
Entrez	1991
Dilution Range	ICC: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.