

Product Datasheet

MCT 1 blocking peptide orb668776

Description	MCT 1 blocking peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	MCT 1
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb668771 (MCT 1 antibody)