

Product Datasheet

CIC-7 blocking peptide orb668777

| | |
|----------------------------|---|
| Description | CIC-7 blocking peptide |
| Conjugation | Unconjugated |
| Tested Applications | HPLC |
| Target | CIC-7 |
| Form/Appearance | Lyophilized |
| Storage | Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles. |
| Note | For research use only. |
| Expiration Date | 12 months from date of receipt. |
| Application Notes | This peptide was used as immunogen to produce orb668772 (CIC-7 antibody) |