

Product Datasheet

Wnt2b Antibody orb654474

Description Wnt2b Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen E.coli-derived human Wnt2b recombinant protein (Position: A7-T391).

Form/Appearance Lyophilized

Concentration Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note For research use only.

Isotype Rabbit IgG

Clonality Polyclonal

MW 60 kDa

Uniprot ID [Q93097](#)

Dilution Range Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date 12 months from date of receipt.

Application Notes Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500 µg/ml.

