

Product Datasheet

HRP Conjugated AffiniPure Donkey Anti-Rabbit IgG(H+L) orb654484

Species/Host	Donkey
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Rabbit IgG (whole molecular)
Form/Appearance	Concentrated, Liquid
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Western Blotting: 0.005-1 µg/ml (ECL detection) Direct ELISA: 0.01-0.2 µg/ml (TMB detection)
Expiration Date	12 months from date of receipt.
Application Notes	Diluted by neutral PBS or TBS. Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.