

Product Datasheet

NUP88 antibody {CF488A} orb221568-CF488A

Description	Rabbit polyclonal antibody to NCR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Equine, Porcine
Conjugation	CF488A
Immunogen	KLH conjugated synthetic peptide derived between 561-640 amino acids of human NUP88
Target	NUP88
Form/Appearance	PBS, 2 mM sodium azide
Concentration	- 100 µg (in 411 &mul;l): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	85 kDa
Uniprot ID	Q99567 (Human); Q8CEC0 (Mouse); O08658 (Rat)
Entrez	4927
Dilution Range	to be determined by end user
Expiration Date	12 months from date of receipt.