

Product Datasheet

FAM3A antibody {CF488} orb11232-CF488

Description	Rabbit polyclonal to FAM3A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF488
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived between 161-220 amino acids of human FAM3A
Target	FAM3A
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	26 kDa
Uniprot ID	P98173 (Human)
Entrez	60343
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.