
Product Datasheet

Rabbit Anti-Acrosin antibody orb500975-CF488A

Description	Rabbit polyclonal antibody to Acrosin
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, FC, WB
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.