

Product Datasheet

Tasosartan orb181331

Description	Tasosartan is a long-acting angiotensin II (AngII) receptor blocker. Its long duration of action has been attributed to its active metabolite enoltasosartan. It is used to treat patients with essential hypertension. For the detailed information of Tasosartan, the solubility of Tasosartan in water, the solubility of Tasosartan in DMSO, the solubility of Tasosartan in PBS buffer, the animal experiment (test) of Tasosartan, the cell experiment (test) of Tasosartan, the in vivo, in vitro and clinical trial test of Tasosartan, the EC50, IC50, and Affinity of Tasosartan, Please contact DC Chemicals..
Note	For research use only.
Formula	C ₂₃ H ₂₁ N ₇ O
MW	411.45914
CAS Number	[145733-36-4]
SMILES	<chem>O=C1N(CC2=CC=C(C3=CC=CC=C3C4=NN=NN4)C=C2)C(N=C5C)=C(C(C)=N5)CC1</chem>
Expiration Date	12 months from date of receipt.