

Product Datasheet

Mouse NLRP3 ELISA Kit (Ready to Use) orb692916

Description	Sandwich Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Mouse
Range	0.156-10ng/mL
Sensitivity	0.055ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8R4B8
NCBI	216799
Entrez	216799
Expiration Date	Please enquire.
Application Notes	serum, plasma, tissue homogenates and other biological fluids.