


Product Datasheet

Vinculin (Loading Control) antibody orb704245

Description	Vinculin (Loading Control) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Human Vinculin
Target	Vinculin (Loading Control)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	124 kDa
Clone ID	5G7
Uniprot ID	P18206
Dilution Range	WB=1:500-2000, ICC=1:50, IF=1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-Vinculin antibody at 1:1000 dilution.