



Product Datasheet

MKLP1 antibody orb704336


Description	MKLP1 antibody
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human MKLP1
Target	MKLP1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	110 kDa
Clone ID	7C9
Uniprot ID	Q02241
Dilution Range	WB=1:500-1000, ICC=1:50, Flow-Cyt=2ug/Test, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of A549 cells with Primary Antibody (green line): Rabbit Anti-MKLP1 antibody; Dilution: 2ul/Test.



Western blot: Lane 1: 293T Cell; Lane 2: A549 Cell; Primary: Anti-MKLP1 at 1:500 dilution.



Western blot: Lane 1: Rat Lymph node; Lane 2: Rat Cerebrum; Lane 3: Rat Testis; Primary: Anti-MKLP1 at 1:1000 dilution.