

Product Datasheet

SPTBN2 antibody orb704489

Description	SPTBN2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Predicted Reactivity	Bovine, Canine, Mouse, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SPTBN2 (71-170/2390 aa)
Target	SPTBN2
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	271 kDa
Uniprot ID	O15020
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.