

Product Datasheet

AAV5 capsid VP2 protein antibody {CF647} orb155552-CF647

Description	Rabbit polyclonal antibody to AAV5 capsid VP2 protein
Conjugation	CF647
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.