

Product Datasheet

Mouse Anti-Chloramphenicol antibody orb156353-CF647

Description	Mouse monoclonal antibody to Chloramphenicol
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA
Storage	Shipped at 4°. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.